



A Turnkey, Fully-Automated, Cloud Computing-Integrated Genomics Solution for Public Health Surveillance of Tomorrow

J. Ng, S. Shokralla, A. Allred and R. Khaksar
Clear Labs, Inc.

INTRODUCTION


- Next generation sequencing (NGS) becoming method of choice for outbreak investigation & antimicrobial resistant (AMR) profiling.¹
- Highly complex & manual NGS workflows typically lead to inaccurate results and costly reruns of tests due to operator-to-operator variability and errors.²
- Automating NGS workflow, from sample extraction to library preparation and sequencing, would significantly improve laboratory efficiency, reduce turnaround times and labor costs.²
- A turnkey, end-to-end, fully automated whole genome sequencing (WGS) solution was developed.
- Cost-effective, gram-agnostic, single touchpoint workflow for surveillance and characterization of bacterial and fungal pathogens.
- Increase productivity and operation efficiency in resource-constrained public health laboratories.^{3,4,5}
- Expedite generation of actionable data during outbreak investigation.^{1,3,4,5}

The Clear Dx™ Automated Platform



METHODS




Overview of the Microbial Surveillance WGS Workflow






¹ Gilchrist, C.A. et al. Whole-genome sequencing in outbreak analysis. *Clin Microbiol Rev.* 2015, 28(3), 541-563.
² Holland, L. & Davies, J.A. Automation in the Life Science Research Laboratory. *Front. Bioeng. Biotechnol.* 2020, 8, 571777.
³ Aboelenen, K., Ivers, P.C. & Abdalrhman, B. Applications of Clear Dx whole genome sequencing system in SARS-CoV-2 diagnostics. *J. Infect. Public Health.* 2022, 15(8), 894-895.
⁴ Lam-Hine, T. et al. Outbreak Associated with SARS-CoV-2 B.1.617.2 (Delta) Variant in an Elementary School – Marin County, California, May–June 2021. *MMWR.* 2021, 70, 1214-1218.
⁵ Ramiah, A. et al. Genomic Surveillance Reveals the Rapid Expansion of the V08 Lineage among Circulating SARS-CoV-2 Omicron Lineages in Southeastern Wisconsin, USA. *Virology.* 2023, 16, 1940.

RESULTS

CDC PulseNet Organisms – Internal Validation



CDC PulseNet Organisms – Field Data




Time and Motion Study

	Clear Dx™	Legacy NGS
Time to results	~27 hrs	4-10 days
Approximate hands-on time	45 min	9-10 hrs
Human touchpoints	1	10+

Other Use Cases: Healthcare-Associated Infection (HAI)

Gram Negative	Gram Positive
<i>Salmonella</i> <i>Escherichia</i> <i>Shigella</i> <i>Vibrio</i> <i>Campylobacter</i> <i>Klebsiella</i> <i>Pseudomonas</i> <i>Acinetobacter</i> <i>Enterobacter</i> <i>Shigella</i> <i>Serratia</i> <i>Citrobacter</i> <i>Raoultella</i> <i>Cronobacter</i>	<i>Providencia</i> <i>Yersinia</i> <i>Proteus</i> <i>Aeromonas</i> <i>Edwardsiella</i> <i>Neisseria</i> <i>Burkholderia</i> <i>Stenotrophomonas</i> <i>Shewanella</i> <i>Kluyvera</i> <i>Pluralibacter</i> <i>Elizabethkingia</i> <i>Eikarella</i>

46 pathogens internally verified using the application to date



CONCLUSION

- Clear Dx™ Microbial Surveillance WGS is a fully automated, end-to-end solution for bacterial and fungal characterization.
- Robust, field-tested and easy-to-use application.
- High quality and reproducible data meeting CDC PulseNet's stringent requirements.
- System accepts dsDNA and amplicons as input material.
- Empower public health laboratories to effectively characterize, monitor and survey pathogens of interest.